


AMRAD Medical



AAU ELITE

Single-detector solution
designed for high-volume facilities


U-ARM RADIOGRAPHIC AND INTEGRATED DR SYSTEM



Large display of arm and detector rotation



Large LCD display with easy touch buttons for accurate positioning



Near-patient controls offer efficient equipment positioning refinements



Position equipment from anywhere in the room



Choose from a variety of optional mobile tables to add to your system.

Peak Performance

The AAU Elite is an advanced, fully-integrated DR System that is driven through worklist commands. Designed to improve the patient experience, the AAU Elite allows the technologists to optimize workflow and improve patient throughput.

Seamless Stitching

Perform LLI and Scoliosis exams with accuracy and ease using the automated equipment and automatic stitching function which comes standard with the AAU Elite.

Smart Investment

Both our AAU Elite and AAU Classic systems offer equal features and earn equal reimbursement to all of the more costly competition . . . all of this at 5-times the warranty. Smart decisions equal smart investments. Be Smart!

**Improve Workflow with
protocol-driven movements.**



Refine equipment positions and remove complicated equipment angles with conveniently located control panel

AMRAD® Medical DR Elite™ Image Capture Workstations



- Simple, clean, and intuitive two-screen workflow: (Work list and Acquire/Review screen)
- Touchscreen or keyboard entries
- Image preview in 2 seconds from exposure
- Integration with most EMR Systems, full MWL capability, MPPS integrations with HIS/RIS providers
- Multiple automated exports that can be user-configured with the single touch of a button
- Complete integration of DAP meters. Readings recorded for export and embedded in DICOM image tags



Features

- Generator Interface Control
- Rejection analysis with user customization of reports
- Embedded Operating System for network safe environments
- Embedded HELP Functions for immediate support services
- Windows 10 Embedded Operating System



Spinal Stitching Study

Standard Advanced Features

- Automatic stitching
- 24/7/365 support with 99% up-time guarantee (1-year included, up to 4 years for purchase)



www.amradmedical.com

SUMMIT INDUSTRIES, LLC
7555 N. Caldwell Avenue, Niles, Illinois 60714
(p) 773-353-4000 | (p) 1-800-729-9729



NSF-ISR

Registered to ISO 13485



LITR041-0723