



SCENARIA View

A powerful premium performance 128-slice CT solution

FUJIFILM
Value from Innovation

Scenaria View CT is designed for radiology departments, EDs, and outpatient centers requiring dependable imaging capabilities, with available advanced modules for interventional CT, extended coverage shuttle scanning for perfusion exams, cardiac CTA, and dual-energy procedures.

Scenaria View features a unique lateral shift table with ± 10 cm shift range and Fujifilm's advanced iterative reconstruction, Intelli-IPV for dose reduction with exceptional low contrast resolution.

Product Highlights



- 80 cm wide bore, 128-Slice CT scanner
- 3-module design fits tight spaces
- New generation low noise HV detector
- Next-generation iterative reconstruction
- Auto lateral shift table
- HiMAR Plus metal artifact reduction
- Rapid Registration and Protocol selection
- SynergyDrive automated workflow shortcuts
- Maximum high-speed scanning 0.35 sec.
- Enhanced Cardiac CTA
- Interventional Needle Guidance
- Dual energy scanning
- Shuttle Scan for brain perfusion exams
- Remote diagnostic maintenance
- Cardio Still Shot intelligent motion suppression coronary imaging
- Integrated Synapse 3D post processing